

Product Information

Date: 2/12/2010

FGS-1000H Test Stand

FGS-1000H Test Stand 1000 lbs (2500N, 500KG) capacity manual test stand.

Designed for convenient testing in both TENSION and COMPRESSION, the FGS-1000H is easy to set-up and operate.

A highly accurate turn-wheel gives a precise 0.1" travel per revolution with a total travel of 6 inches. In addition, it's rigged and rugged construction provides repeatable and reliable testing, time after time. The FGS-1000H features fully adjustable height ranges in 13 one inch increments accommodating many different sample heights.

With counter sunk (force gauge) mounting holes, a self aligning, pre-tapped sliding sample base plate, the FGS-1000H provides a welcome alternative to utilizing existing force gauges, grips attachments and accessories.

